

Product Data Sheet

Recombinant Human EDIL3 Protein

Catalog #	Source	Reactivity	Applications
CST2375	E. coli	H	
Description	Recombinant Human EDIL3 Protein is produced by E. coli expression system and the target gene encoding Glu321-Arg470 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	EDIL3		
Alternative Names	DEL1; EGF-like repeat and discoidin I-like domain-containing protein 3; Developmentally-regulated endothelial cell locus 1 protein; Integrin-binding protein DEL1		
Entrez Gene	10085		
SwissProt	O43854		
Purity			
Chemical Structure	EPLGMKSGHI QDYQITASSI FRTLNMMDMFT WEPRKARLDK QGKVNAWTSG HNDQSQWLQV DLLVPTKVTG IITQGAKDFG HVQFVGSYKL AYSNDGEHWT VYQDEKQRKD KVFQGNFDND THRKNVIDPP IYARHIRILP WSWYGRITLR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com