

Product Data Sheet

Recombinant Human FLRT2 Protein

Catalog #	Source	Reactivity	Applications
CST2376	E. coli	H	
Description	Recombinant Human FLRT2 Protein is produced by E. coli expression system and the target gene encoding Arg41-Thr230 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FLRT2		
Alternative Names	KIAA0405; Leucine-rich repeat transmembrane protein FLRT2; Fibronectin-like domain-containing leucine-rich transmembrane protein 2		
Entrez Gene	23768		
SwissProt	O43155		
Purity			
Chemical Structure	RCDRNFVYCN ERSLTSVPLG IPEGVTVLVL HNNQINNAGF PAELHNVQSV HTVYLYGNQL DEFPMNLPKN VRVLHLQENN IQTISRALA QLLKLEELHL DDNSISTVGV EDGAFREAIS LKLLFLSKNH LSSVPVGLPV DLQELRV DEN RIAVISDMAF QNLTSLERLI VDGNLLTNKG IAEGTFSHLT		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com