

Product Data Sheet

L-Amino Acid Oxidase

Catalog #	Source	Reactivity	Applications
CZM1172	Microorganism		
Description	Recombinant L-Amino Acid Oxidase from Microorganism		
Specificity	>4 U/mg Protein, >2 U/mg powder		
Form	Light yellow lyophilized powder		
Alternative Names	EC 1.4.3.2; laao; il4i1; l-amino-acid oxidase; l-aao; escapin; head kidney and gill; dolabellatin; l-phenylalanine oxidase; akbu-laao; m-lao		
CAS Number	9000-89-9		
Molecular Weight			
Purity	>90%		
Unit Definition	One unit is defined as the amount of enzyme will oxidize 1 micromole L-leucine per minute at pH 7.6, 25°C.		
Storage/Stability	Shipped at 4°C. Store at -20°C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com