

Product Data Sheet

endo-1,3-beta-glucanase

Catalog #	Source	Reactivity	Applications
CZM1176	Microorganism		
Description	Recombinant endo-1,3-beta-glucanase from Microorganism		
Specificity	>10 U/mg Protein, >1 U/mg powder		
Form	White lyophilized powder		
Alternative Names	EC 3.2.1.39; glucan endo-1,3-beta-D-glucosidase; (1->3)-beta-glucan endohydrolase; laminarinase		
CAS Number	9044-93-3		
Molecular Weight			
Purity	>90%		
Unit Definition	One unit is defined as the amount of enzyme will liberate 1 micromole glucose from laminarin per minute at pH 5.0, 37°C.		
Storage/Stability	Shipped at 4°C. Store at -20°C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com